# Mole®

## GeneMole® Manuel d'utilisation

MGM-102-005

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## 13. GeneMole® Troubleshooting Guide

| GeneMole® Troubleshooting Guide |  |          |
|---------------------------------|--|----------|
| Issue                           | Problem  | Refer to |
| 1                               | The instrument seems dead (no pilot light, nothing on screen).   | Table 1  |
| 2                               | The instrument pilot light is on, but there is nothing on the screen   | Table 2  |
| 3                               | The screen shows the Windows CE desktop.   | Table 3  |
| 4                               | The screen picture is normal for GeneMole®, but the screen does not react to touch (frozen screen).  | Table 4  |
| 5                               | The instrument will not start a run, error message “Door opened”.  | Table 5  |
| 6                               | The instrument stops with the error message “Crash sensor activated” or “The pipette head has encountered an unexpected obstacle”.   | Table 6  |
| 7                               | The instrument displays the error message “Y1 limit switch activated”.   | Table 7  |
| 8                               | The instrument displays the error message “X1 limit switch activated”.   | Table 8  |
| 9                               | The instrument displays the error message “X2 limit switch activated”.   | Table 9  |
| 10                              | The instrument displays the error message “Y2 limit switch activated”.   | Table 10 |
| 11                              | The instrument displays the error message “P limit switch activated”.  | Table 11 |
| 12                              | The pipette head loses tips.   | Table 12 |
| 13                              | The pipettes seem to leak.   | Table 13 |
| 14                              | When volume calibration is performed, the new data are not accepted because no calibration file exists for the present pipette head.                                       | Table 14 |
| 15                              | When volume calibration is performed, the new data are not accepted because the variation is too great or because the new data deviate too much from the factory settings. | Table 15 |
| 16                              | One axis does not move, although the instrument does not give an error message.  | Table 16 |
| 17                              | Liquid overflow/spills.  | Table 17 |
| 18                              | Instrument misses positions, does not pick up tips or lifts items from the work tray.  | Table 18 |
| 19                              | Variable elution volume.   | Table 19 |
| 20                              | Reduced yield and/or quality.  | Table 20 |

**Table 1**

| <b>Problem</b>  | <b>Cause / Remedy</b>   |
|---|---|
| <b>The instrument seems dead (no pilot light, nothing on screen).</b> | <ol style="list-style-type: none"><li>1. Check that the instrument is turned on, if not turn on.</li><li>2. Check that the instrument is connected to the mains, if not, connect.</li><li>3. Check that the mains is in working order (use a lamp or other appliance), if not, switch to another mains socket.</li><li>4. Check the instrument fuse. If blown, replace (refer to User Manual section 3).</li><li>5. If all of the above is OK, turn the instrument off, leave it for 5 minutes and turn it on again. If the instrument does not come alive within three minutes after being turned on this time, call for service giving the reference A01.</li></ol> |

**Table 2**

| <b>Problem</b>   | <b>Cause / Remedy</b>   |
|--|---|
| <b>The instrument pilot light is on, but there is nothing on the screen.</b> | <ol style="list-style-type: none"><li>1. Turn the instrument off, then on again.<ol style="list-style-type: none"><li>a) If the instrument starts directly, go ahead and use as normal.</li><li>b) If the instrument shows the splash screen for longer than normal, and states that the system is being updated, recalibration may be needed before use. Do a test run and if there seems to be a problem, call for service giving the reference A02.</li></ol></li><li>2. If the instrument does not restart, turn the instrument off, leave it for 5 minutes and turn it on again.<ol style="list-style-type: none"><li>a) If the instrument starts directly, go ahead and use as normal.</li><li>b) If the instrument does not come alive within three minutes after being turned on this time, call for service giving the reference A03.</li><li>c) If the instrument shows the splash screen for longer than normal, and states that the system is being updated, recalibration may be needed before use. Do a test run and if there seems to be a problem, call for service giving the reference A02.</li></ol></li></ol> |

**Table 3**

| Problem  | Cause / Remedy   |
|--|--|
| <p><b>The screen shows the Windows CE desktop.</b></p> | <ol style="list-style-type: none"> <li>1. Turn the instrument off, then directly on again.               <ol style="list-style-type: none"> <li>a) If the instrument starts directly, go ahead and use as normal.</li> <li>b) If the instrument shows the splash screen for longer than normal, and states that the system is being updated, recalibration may be needed before use. Do a test run and if there seems to be a problem, call for service giving the reference A02.</li> </ol> </li> <li>2. If the instrument does not restart, turn the instrument off, leave it for 5 minutes and turn it on again.               <ol style="list-style-type: none"> <li>a) If the instrument starts directly, go ahead and use as normal.</li> <li>b) If the instrument does not come alive within three minutes after being turned on this time, call for service giving the reference A03.</li> <li>c) If the instrument shows the splash screen for longer than normal, and states that the system is being updated, recalibration may be needed before use. Do a test run and if there seems to be a problem, call for service giving the reference A02.</li> </ol> </li> </ol> |

**Table 4**

| Problem   | Cause / Remedy  |
|---|---|
| <p><b>The screen picture is normal for GeneMole®, but the screen does not react to touch (frozen screen).</b></p> | <ol style="list-style-type: none"> <li>1. Turn the instrument off, then directly on again.               <ol style="list-style-type: none"> <li>a) If the instrument starts directly, go ahead and use as normal.</li> <li>b) If the instrument shows the splash screen for longer than normal, and states that the system is being updated, recalibration may be needed before use. Do a test run and if there seems to be a problem, call for service giving the reference A02.</li> </ol> </li> <li>2. If the instrument does not restart, turn the instrument off, leave it for 5 minutes and turn it on again.               <ol style="list-style-type: none"> <li>a) If the instrument starts directly, go ahead and use as normal.</li> <li>b) If the instrument does not come alive within three minutes after being turned on this time, call for service giving the reference A03.</li> </ol> <p>If the instrument shows the splash screen for longer than normal, and states that the system is being updated, recalibration may be needed before use. Do a test run and if there seems to be a problem, call for service giving the reference A02.</p> </li> </ol> |

Table 5

| Problem   | Cause / Remedy  |
|---|---|
| <p><b>The instrument will not start a run, error message “Door opened”.</b></p> | <ol style="list-style-type: none"> <li>1. Check that the door is properly closed, if not, close it.</li> <li>2. If the door is properly closed and the error message still appears, try to depress the door switches by hand. If the error message disappears, the problem is a maladjusted door. This information should be given when calling for service giving the reference A04.</li> <li>3. If depressing the door switches by hand does not remove the error message “Door opened”, the error is in the electronics system. This information should be given when calling for service giving the reference A05.</li> </ol> |

Table 6

| Problem  | Cause / Remedy  |
|--|---|
| <p><b>The instrument stops with the error message “Crash sensor activated” or “The pipette head has encountered an unexpected obstacle”.</b></p> <p>This means that the pipette head has encountered more resistance than normal when moving down.</p> | <ol style="list-style-type: none"> <li>1. Check for the presence of obstacles that may have stopped the pipette head from moving down. If found, remove them and restart the instrument.</li> <li>2. Make sure that only plastic ware recommended by Mole Genetics is used. Other brands of plastic, though superficially similar, may have different dimension and cause the instrument to crash.</li> <li>3. Make sure that sample tubes have been pressed all the way into their holes (they should rest on the small collar). If this is not the case, press them down.</li> <li>4. Make sure that the black plate is in position between the <b>MoleStrips™</b> and the strip holder. If the <b>MoleStrips™</b> are placed directly into the strip holder, they will sit lower than normal, and the instrument will expect to find the bottom of the sample tubes correspondingly lower.</li> <li>5. If no obvious obstacles are found, restart the instrument. If the problem persists, call for service giving the reference A06.</li> </ol> |

**Table 7**

| Problem  | Cause / Remedy   |
|--|--|
| <p><b>The instrument displays the error message “Y1 limit switch activated”.</b></p> <p>The Y1 limit switch is related to the vertical movement of the pipette head and this problem is similar to “Crash sensor activated”. At some steps resistance to movement is expected, and the crash sensor is turned off. If the resistance at those steps is even greater than expected, the y1 motor may lose its position, causing it to go too far when it returns to its upper position.</p> | <ol style="list-style-type: none"> <li>1. Check for the presence of obstacles that may have stopped the pipette head from moving down. If found, remove them and restart the protocol.</li> <li>2. Check that the <b>MoleStrips™</b> have been properly seated in the strip holder (with black plate). If the <b>MoleStrips™</b> are placed higher at one end than the other, this may cause this error.</li> <li>3. If no obvious obstacles are found, restart the instrument. If the problem persists, call for service giving the reference A07.</li> </ol> |

**Table 8**

| Problem  | Cause / Remedy  |
|--|---|
| <p><b>The instrument displays the error message “X1 limit switch activated”.</b></p> <p>The X1 limit switch is related to the horizontal movement of the pipette head.</p> | <ol style="list-style-type: none"> <li>1. Check for obstacles to sideways movement of pipette head. If found, remove them.</li> <li>2. Check if the pipette head has been unable to detach pipette tips at the correct stage, since this may cause the pipette head to collide with the magnet arm. Possible causes for non detachment of tips is that tips from previous runs have been left in the waste container or that something has fallen into the waste container.</li> <li>3. If no obvious obstacles are found, restart the instrument. If the problem persists, call for service giving the reference A08.</li> </ol> |

**Table 9**

| Problem  | Cause / Remedy   |
|--|--|
| <p><b>The instrument displays the error message “X2 limit switch activated”.</b></p> <p>The X2 limit switch is related to the horizontal movement of the magnet arm.</p> | <ol style="list-style-type: none"> <li>1. Check for obstacles to sideways movement of magnet arm. If found, remove them.</li> <li>2. Check if the pipette has been unable to detach pipette tips at the correct stage, since this may cause the pipette head to collide with the magnet arm. Possible causes for non detachment of tips is that tips from previous runs have been left in the waste container or that something has fallen into the waste container.</li> <li>3. If no obvious obstacles are found, restart the instrument. If the problem persists, call for service giving the reference A09.</li> </ol> |

**Table 10**

| Problem  | Cause / Remedy  |
|--|---|
| <p><b>The instrument displays the error message “Y2 limit switch activated”.</b></p> <p>The Y2 limit switch is related to the vertical movement of the magnet arm.</p> | <ol style="list-style-type: none"><li>1. Check for the presence of obstacles that may have stopped the magnet arm from moving down. If found, remove them and restart the instrument.</li><li>2. If no obvious obstacles are found, restart the instrument. If the problem persists, call for service giving the reference A10.</li></ol> |

**Table 11**

| Problem   | Cause / Remedy  |
|---|---|
| <p><b>The instrument displays the error message “P limit switch activated”.</b></p> | <ol style="list-style-type: none"><li>1. Check if the sample has been viscous enough to completely block the pipette tips. If this is the case, use less sample and/or avoid clots.</li><li>2. Check that nothing blocks the nozzles of the pipette head.</li><li>3. If no obvious obstacles are found, restart the instrument. If the problem persists, call for service giving the reference A11.</li></ol> |

**Table 12**

| Problem                                    | Cause / Remedy  |
|--|---|
| <p><b>The pipette head loses tips.</b></p> | <ol style="list-style-type: none"><li>1. Call for service giving the reference A12.</li></ol> |

**Table 13**

| Problem  | Cause / Remedy   |
|--|--|
| <p><b>The pipettes seem to leak.</b></p> <p>This may be a real leakage due to a damaged pipette head, but there are also several steps in the protocols where the tips have been deep in viscous liquid. Such liquid may then run down the outside of the tips and form a drop at the end of the tip. This does not in any way compromise the run.</p> | <ol style="list-style-type: none"> <li>1. Run the protocol “Leakage test” as described in section 10.2.2 in the User Manual to check if the tips are leaking.               <ol style="list-style-type: none"> <li>a) If no drops have formed after 5 minutes, leakage is not a problem.</li> <li>b) If drops have formed at one or more tips after 5 minutes call for service, indicating on how many tips drops have formed giving the reference A13.</li> </ol> </li> </ol> |

**Table 14**

| Problem  | Cause / Remedy  |
|--|---|
| <p><b>When volume calibration is performed, the new data are not accepted because no calibration file exists for the present pipette head.</b></p> | <ol style="list-style-type: none"> <li>1. Open the pipette head cover and check that the switch plate depresses the two pipette head recognition switches.               <ol style="list-style-type: none"> <li>a) If the switch plate of the pipette head does not press down the two switches, twist it until it does. This may require loosening the pipette head holding nut.</li> <li>b) If the switch plate of the pipette head does press down the two switches, call for service giving the reference A14.</li> </ol> </li> </ol> |

**Table 15**

| Problem  | Cause / Remedy  |
|--|---|
| <p><b>When volume calibration is performed, the new data are not accepted because the variation is too great or because the new data deviate too much from the factory settings.</b></p> | <ol style="list-style-type: none"> <li>1. The most common reason for a failed volume calibration is typing error. Restart the instrument and do the calibration again. Note that if the instrument is not restarted after a failed calibration, it will continue to fail, even if the values are within the predefined limits               <ol style="list-style-type: none"> <li>a) If the instrument accepts new calibration values after restart, the problem is solved.</li> <li>b) If the instrument continues to reject the new calibration values, it is probable that the pipette head has to be replaced, call for service giving the reference A15.</li> </ol> </li> </ol> |

Table 16

| Problem   | Cause / Remedy  |
|---|---|
| One axis does not move, although the instrument does not give an error message. | 1. This is always a matter for a service engineer, call for service giving the reference A16. |

Table 17

| Problem                | Cause / Remedy  |
|------------------------|---|
| Liquid overflow/spills | <ol style="list-style-type: none"> <li>1. Check that the right sample volume was applied, including whether the pipette used to dispense the sample is correctly calibrated.</li> <li>2. Check that there is nothing in the samples that may clog the pipette tips (e.g. blood clots, undigested tissue remains, hair).</li> <li>3. Make sure that only plastic ware recommended by Mole Genetics is used. Other brands of plastic, though superficially similar, may have different dimension and cause liquid overflow.</li> <li>4. Make sure that sample tubes have been pressed all the way into their holes (they should rest on the small collar). If this is not the case, press them down.</li> <li>5. Make sure that the black plate is in position between the <b>MoleStrips™</b> and the strip holder. If the <b>MoleStrips™</b> are placed directly into the strip holder, they will sit lower than normal, and the instrument will expect to find the bottom of the sample tubes correspondingly lower. This may cause liquid overflow.</li> <li>6. Run the protocol “Leakage test” as described in section 10.2.2 in the User Manual to check if the tips are leaking. <ol style="list-style-type: none"> <li>a) If no drops have formed after 5 minutes, leakage is not a problem</li> <li>b) If drops have formed at one or more tips after 5 minutes call for service, indicating on how many tips drops have formed giving the reference A13</li> </ol> </li> <li>7. If the problem persists, even though there is no leakage from the pipette, call for service giving the reference A17.</li> </ol> |

Table 18

| Problem   | Cause / Remedy   |
|---|--|
| <p><b>Instrument misses positions, does not pick up tips or lifts items from the work tray.</b></p> | <ol style="list-style-type: none"> <li>1. Make sure the tray is correctly situated. The two guide pins at the front of the instrument floor should go through the corresponding holes in the tray. Correct if necessary.</li> <li>2. If the problem persists, it is an indication that a recalibration of the factory set positions is needed, call for service giving the reference A18.</li> </ol> |

Table 19

| Problem                                | Cause / Remedy   |
|--|--|
| <p><b>Variable elution volume.</b></p> | <ol style="list-style-type: none"> <li>1. Check that the <b>MoleStrips™</b> have been properly seated in the strip holder (with black plate). If <b>MoleStrips™</b> are higher at one end than the other, this may cause this error.</li> <li>2. Check that there is liquid left in the well where the elution buffer was aspirated from. It is necessary to use a pipette to check this, the profile of the lower part of the well makes visual inspection difficult.             <ol style="list-style-type: none"> <li>a) If there is no elution buffer left, visually inspect some of the other <b>MoleStrips™</b> from the same kit. If elution buffer volumes seem uneven, this is a likely cause.</li> <li>b) If the <b>MoleStrips™</b> seem in order, perform a volume calibration.</li> </ol> </li> <li>3. If there is elution buffer left in the strip, run the protocol “Leakage test” to check if the tips are leaking.             <ol style="list-style-type: none"> <li>a) If no drops have formed after 5 minutes, leakage is not a problem.</li> <li>b) If drops have formed at one or more tips after 5 minutes call for service, indicating on how many tips drops have formed giving the reference A13.</li> </ol> </li> <li>4. If the problem persists even though there is no leakage, no problem with the kits and a correct volume calibration has been performed, call for service giving the reference A19.</li> </ol> |

**Table 20**

| Problem                             | Cause / Remedy  |
|-------------------------------------|---|
| <b>Reduced yield and/or quality</b> | <ol style="list-style-type: none"><li data-bbox="715 376 1469 439">1. If the problem seems related to any of the problems described above, try out the solutions described for that problem.</li><li data-bbox="715 456 1501 519">2. If the problem is unrelated to any problem described above, refer to the documentation of the kit used.</li><li data-bbox="715 537 1501 631">3. If the problem is unrelated to any problem described above and neither seem to be related to the kit, or if it affects more than one type of kit, call for service giving the reference A20.</li></ol> |

## Notes

## Notes

## Notes

Mole Genetics AS  
Vollsveien 13D, N-1366 Lysaker, Norvège  
Phone: + 47 6710 1820, Fax: + 47 6710 1821

[www.molegenetics.com](http://www.molegenetics.com)